

Cable and Wireless Fitting Guide

Fitting with ReSound Aventa® 3									
 17292002 (Airlink 2, GN Hearing) 17292000 (Airlink, GN Hearing)	 9022 907 69019 (HI-PRO CABLE/RIGHT/CS44) 9022 907 69029 (HI-PRO CABLE/LEFT/CS44) 3122 150 73781 (FLEX CONNECTOR CLIP CS63) 19181700 (STRIP,FLEX,CS63,PRG)	 15206700 (Programming cable CS63, Right) 15206701 (Programming cable CS63, Left) 19181700 (STRIP,FLEX,CS63,PRG)	 9022 907 69019 (HI-PRO CABLE/RIGHT/CS44) 9022 907 69029 (HI-PRO CABLE/LEFT/CS44) 19533900 (ADAPTER,FLEX,3 PIN,CS53)	 9022 907 69019 (HI-PRO CABLE/RIGHT/CS44) 9022 907 69029 (HI-PRO CABLE/LEFT/CS44) 18041400 (PROGRAMMING ADAPTOR,8PIN)	 9022 907 69019 (HI-PRO CABLE/RIGHT/CS44) 9022 907 69029 (HI-PRO CABLE/LEFT/CS44) 17598200 (PRGM.SHOE ASM,W/ CS44,DELHI)	 18523700 (IIC 10A PROG PILL CABLE, LEFT) 18523800 (IIC 10A PROG PILL CABLE, RIGHT) 17167800 (PROGRAMMING PILL)	 9022 907 69019 (HI-PRO CABLE/RIGHT/CS44) 9022 907 69029 (HI-PRO CABLE/LEFT/CS44) 17167800 (PROGRAMMING PILL)	 9022 907 69019 (HI-PRO CABLE/RIGHT/CS44) 9022 907 69029 (HI-PRO CABLE/LEFT/CS44) 19215000 (KIT,PROG,SUPER POWER, LO90)	
ReSound LiNX ² ™									
RIE WL - LSx62	•				•				
RIE WL - LSx61	•	•	•	•	•				
BTE WL - LSx77, LSx88	•				•				
Custom WL	•	•	•	•	•				
Custom non-WL		•	•	•	•				
Custom IIC							•		
ReSound ENZO™									
BTE - EOx98	•								•
ReSound LiNX™									
RIE WL - LNx61	•	•	•	•	•				
BTE WL - LNx77, LNx88	•				•				
ReSound Verso™ & ReSound Alera®									
RIE WL - VOx62 / ALx62	•				•				
RIE WL - VOx61 / ALx61	•	•	•	•	•				
RIE - VOx60 / ALx60		•	•	•	•				
BTE WL - VOx67, VOx77, VOx88	•				•				
BTE WL - ALx67, ALx77/87	•				•				
Custom WL	•	•	•	•	•				
Custom		•	•	•	•				
Custom IIC - VOx10-C							•		
ReSound Lex™									
RM - LXx00-M								•	
ReSound Enya™									
BTE WL - EYx67, EYx77, EYx88	•				•				
RIE WL - EYx62	•				•				
Custom WL	•	•	•	•	•				
Custom		•	•	•	•				
ReSound Magna™ & ReSound Vea™									
BTE - MGx90						•			
BTE - VE377					•				
RIE - VE362					•				
BTE - VE360, VEx70, VEx80						•			
Custom		•	•	•	•				

Please note that all ReSound hearing instruments for both wireless and non-wireless **fittings require a fresh battery when fitted with Aventa 3.9.**